powered by

NEURAC<sup>®</sup>

# REDCORD MINI

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Opplever du problemer med monteringen, kontakt en forhandler.

Sikkerhetstest av Redcord Mini

## Read and understand the following information before using Redcord Mini.

Redcord Mini is a portable and functional rehabilitation and training device. Intended use of the Redcord Mini is for rehabilitation, injury prevention and for functional strength training. The Redcord Mini enables a fullbody work-out whether you are at home, outdoors or traveling.

Failure to follow these guidelines may increase risk for injury and or reduce lifetime expectancy of the Mini device (Normal use = 10 years). Detailed information on [www.redcord.com](http://www.redcord.com)

- Consult your Physician before using the equipment
  - Inspect the equipment before use
  - Do not use damaged equipment, replace worn out equipment immediately
  - Make sure installation has been done correctly
  - Perform the safety test before first use and then on a monthly basis by putting one foot in the equipment and bounce carefully with full body weight.
  - Do not stand on the rope when exercising. This can result in the ropes sliding out of the locking mechanism and you falling down
- Redcord takes no responsibility for any personal or physical damages that may arise due to incorrect installation and/or incorrect use of Redcord equipment.
  - Redcord is not responsible for any personal or physical damage that may arise due to the installation or due to the ceiling
  - Do not put the rope around your neck

- Follow instructions from the health care professional when performing exercise
  - Stop exercising when experiencing pain or discomfort
  - Be aware that sweat increases the chance of slippery equipment and floor
  - Only use the equipment for its intended purpose
  - Respect the weight load limit, max load marked on the device, 100kg (220 lb) per rope
  - Avoid chock loads - chock loads may damage your device permanently
  - Redcord Mini shall not be used as a swing - repetitive movements will involve much more wear and tear than any normal use and may reduce the strength permanently
  - Children using the equipment should always be supervised
- Store in a cool and dry place with good air circulation when not in use
  - The device can be washed with a soft brush and lukewarm soapy water
  - Washed device shall not be dried in sunlight or in a dryer
  - Note that prolonged exposure to sunlight, dampness, bleach, extreme heat or cold and wrong use can weaken the strength of the device
  - Inspect the ceiling brackets by grabbing the white ceiling brackets and try to move the brackets in vertical and horizontal direction
  - The ceiling brackets should be stable in all directions, tighten the screws if necessary.
  - Document your safety test of the device

## Installation guide

Seek professional advice about the condition of the ceiling prior to the installation of the device. Redcord (nor its distributors) takes no responsibility for the suitability or condition of the ceiling in which the device is installed.

### Ceiling heights

Recommended installation height for Redcord Mini is between 220cm (7'2") and 260cm (9'6").

### Distance between the ropes

The recommended distance between the ropes is 60 - 90 cm (24-36in)

### Ceiling suspension to Redcord Mini

### Wooden ceilings

Make sure you are installing the screws in the ceiling joists and not just in panels or wallboard. Use the accompanying installation screws and pre-drill with a 6mm bit. Test the installation using the safety test.

### Concrete ceilings

Make sure that the ceiling is concrete and not made of a porous material with a poor load capacity, (Siporex, Ytong, Leca). Drill two holes 70mm deep with a 12mm masonry bit. Ensure that the hole is deep enough for the plastic plug to fill the entire hole. Screw in the installation screws and test by performing the safety test.

### Porous ceilings

Do not install Redcord Mini in lightweight aerated concrete (Ytong, Siporex or Leca) ceilings

### Other installation options

Redcord Mini can be used anywhere you can locate a supporting beam of sufficient height. Make sure the beam will support your fully bodyweight by always performing the safety test before performing the exercises

## Installation/Disassembly

Thread the rope through the provided hooks or over a tree branch, bar or beam. Thread the end of the rope through the locking mechanism from above. Test the load tolerance by performing the safety test.

To pull out the adjustable rope from the locking mechanism, push the grey handle while pulling the rope out. This will prevent the rope to get locked inside the locking mechanism. The grey handle can be released when the rope is fully out of the locking mechanism

## Adjusting the rope length

Hold the rope with the strap in one hand. Pull the rope towards you and adjust up or down until you reach the desired height. When hanging parallel, the ropes will lock automatically upon load.

**WARNING:** Do not stand on the rope when exercising with Redcord Mini. This can result in the ropes sliding out of the locking mechanism and you falling forward.

## Safety test

**After installing the ceiling suspension and before you start exercising with Redcord Mini, you must test the installation and the load capacity**

Attached the device in the suspension. Lower the strap to about 15 cm (6 inches) above the floor level. Place one foot in one strap and try to bounce up and down using your full body weight. Repeat the test on the other ceiling suspension. You have now tested that each anchor point tolerates a load far exceeding your body weight. The ceiling brackets are ready to use for training. If you have problem with installation, contact your distributor.

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